

## Rapid Fentanyl Urine Test

For *in vitro* diagnostic use only.

**CLIA WAIVED**

Medical and other professional *in vitro* diagnostic use labeling.

### INTENDED USE

The *EKLA Rapid Fentanyl Urine Test* is a rapid visual immunoassay for the qualitative, presumptive detection of Fentanyl in human urine specimens at a cut-off concentration of 1.0 ng/mL. It is for *in vitro* diagnostic use only.

The tests provide only preliminary test results. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. GC/MS or LC/MS is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary result is positive.

### Contents of the kit

#### Materials Provided

- Test devices
- Disposable pipettes
- Package insert

#### Materials Required but Not provided

- Timer
- Specimen collection container

### PRECAUTIONS

- For *in vitro* diagnostic use only.
- Do not use after the expiration date indicated on the package.
- Do not use the test if the foil pouch or canister is damaged.
- Do not reuse tests.
- Avoid cross-contamination of specimens by using a new specimen collection container for each specimen obtained.
- Read the entire procedure carefully prior to testing.
- Humidity and temperature can adversely affect results.

### STORAGE AND STABILITY

- The kit should be stored at 36-86°F (2-30°C) until the expiry date printed on the sealed pouch.
- The test must remain in the sealed pouch or closed canister until use.
- Do not freeze.
- Kits should be kept out of direct sunlight.

### PROCEDURE


Bring tests to room temperature (59-86°F (15-30°C)) before use.

#### 1. Collect Urine Sample

Collect urine in a clean and dry container.



2. Remove the test from its sealed pouch, and place it on a clean, level surface.
3. Using the provided disposable pipette, transfer 3 drops of specimen to the specimen well (S) of the device and start the timer.  
**Do not add any solution to the result area.**  
As the test begins to work, color will migrate across the membrane.
4. Wait for the colored band(s) to appear. The result should be read at 5 minutes. Do not interpret the result after 10 minutes.



### INTERPRETATION OF RESULTS



**POSITIVE:** Only one colored band appears, in the control region (C). No apparent colored band appears in the test region (T).



**NEGATIVE:** Two colored bands appear on the membrane. One band appears in the control region (C) and another band appears in the test region (T).



**INVALID:** Control band fails to appear. Results from any test which has not produced a control band at the specified read time must be discarded. Please review the procedure and repeat with a new test. If the problem persists, discontinue using the kit immediately and contact your local distributor.

#### NOTE:

1. The intensity of color in the test region (T) may vary depending on the concentration of analytes present in the specimen. Therefore, any shade of color in the test region should be considered negative. Note that this is a qualitative test only, and cannot determine the concentration of analytes in the specimen.
2. Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control band failure.

### LIMITATIONS

1. The *EKLA Rapid Fentanyl Urine Test* should be only used for the qualitative detection of Fentanyl.
2. This assay provides a preliminary analytical test result only.
3. There is a possibility that technical or procedural errors as well as other substances and factors may interfere with the test and cause false results.
4. Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. Therefore, please preclude the possibility of urine adulteration prior to testing.
5. A positive result indicates the presence of a Fentanyl only, and does not indicate or measure intoxication.
6. A negative result does not at any time rule out the presence of Fentanyl in urine, as they may be present below the minimum detection level of the test.
7. This test does not distinguish between Fentanyl and certain medications.

### QUALITY CONTROL

1. Internal procedural controls are included in the test. A colored band appearing in the control region (C) is considered an internal procedural control, confirming sufficient specimen volume and correct procedural technique.
2. External controls are not supplied with this kit. It is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.
3. User should follow federal, state and local guidelines for testing quality control materials. Laboratories should comply with all federal state, and local laws, as well as all guidelines and regulations.

### PERFORMANCE CHARACTERISTICS

#### A. Accuracy

The accuracy of the test was confirmed by testing 82 clinical urine specimens in parallel with LC-MS. The results are shown below.

Test	Negative	Low Negative by LC/MS (less than -50%)	Near Cutoff Negative by LC/MS (Between -50% and cutoff)	Near Cutoff Positive by LC/MS (Between cutoff and +50%)	High Positive by LC/MS (greater than +50%)	
						Operator 1
Operator 1	Positive	0	0	2	16	22
	Negative	5	16	19	2	0
Operator 2	Positive	0	0	2	16	22
	Negative	5	16	19	2	0
Operator 3	Positive	0	0	3	17	22
	Negative	5	16	18	1	0

% agreement among positives is 95.8%

% agreement among negatives is 94.4%

#### B. Specificity and Cross-reactivity

The minimum concentration required to obtain a positive result of the drug metabolites and other components for the test is listed below.

Compounds	Minimum concentration required to obtain a positive result (ng/mL)	% Cross-Reactivity
Acetyl fentanyl	5	20.00
Acrylfentanyl	5	20.00
ω-1-Hydroxyfentanyl	50000	0.002
Isobutyryl fentanyl	5	20.00
Ocfentanil	100	1.00
Butyryl fentanyl	25	4.00
Furanyl fentanyl	10	10.00
Valeryl fentanyl	50	2.00
(±) β-hydroxythiofentanyl	5	20.00
4-Fluoro-isobutyrylfentanyl	50	2.00
Para-fluorobutyryl fentanyl	25	4.00

Para-fluoro fentanyl	25	4.00
(±)-3-cis-methyl fentanyl	50	2.00
Carfentanil	10000	0.01
Despropionyl fentanyl (4-ANPP)	50000	0.002
Sufentanil	>10000	<0.01
Alfentanil	>100000	<0.001
Remifentanil	>10000	<0.01
Norfentanyl	>100000	<0.001
Acetyl norfentanyl	>100000	<0.001
Norcarrfentanil	>10000	<0.01

The opioids compounds at the concentration listed below are verified to be no cross-reactivity for the test.

Compounds	Concentration
6-Acetyl morphine	100 ug/mL
Amphetamine	100 ug/mL
Buprenorphine	100 ug/mL
Buprenorphineglucuronide	100 ug/mL
Codeine	100 ug/mL
Dextromethorphan	100 ug/mL
Dihydrocodeine	100 ug/mL
EDDP	100 ug/mL
EMDP	100 ug/mL
Fluoxetine	100 ug/mL
Heroin	100 ug/mL
Hydrocodone	100 ug/mL
Hydromorphone	100 ug/mL
Ketamine	100 ug/mL
Levorphanol	100 ug/mL
Meperidine	100 ug/mL
Methadone	100 ug/mL
Morphine	100 ug/mL
Morphine-3-glucuronide	100 ug/mL
Naloxone	100 ug/mL
Naltrexone	100 ug/mL
Norbuprenorphine	100 ug/mL
Norcodeine	100 ug/mL

Norketamine	100 ug/mL
Normeperidine	100 ug/mL
Normorphine	100 ug/mL
Noroxycodone	100 ug/mL
Oxycodone	100 ug/mL
Oxymorphone	100 ug/mL
Pentazocine (Talwin)	100 ug/mL
Pipamperone	100 ug/mL
Risperidone	100 ug/mL
Tapentadol	100 ug/mL
Thioridazine	100 ug/mL
Tilidine	100 ug/mL
Tramadol	100 ug/mL
Tramadol-O-Desmethyl	100 ug/mL
Tramadol-N-Desmethyl	100 ug/mL
Trazodone	100 ug/mL

### C. Precision

Result Drug Conc. % of cutoff	Lot 1	Lot 2	Lot 3
-100%	60(-)/ 0(+)	60(-)/ 0(+)	60(-)/ 0(+)
-75%	60(-)/ 0(+)	60(-)/ 0(+)	60(-)/ 0(+)
-50%	60(-)/ 0(+)	60(-)/ 0(+)	60(-)/ 0(+)
-25%	58(-)/2(+)	58(-)/2(+)	57(-)/ 3(+)
cutoff	35(+)/ 25(-)	36(+)/24(-)	35(+)/ 25(-)
+25%	60(+)/0(-)	60(+)/0(-)	60(+)/0(-)
+50%	60(+)/0(-)	60(+)/0(-)	60(+)/0(-)
+75%	60(+)/0(-)	60(+)/0(-)	60(+)/0(-)
+100%	60(+)/0(-)	60(+)/0(-)	60(+)/0(-)

### D. Effect of Urine Density and pH

Eight urine samples of different urine density were tested. The results showed that varying urine density did not affect the test results.

Six urine samples of different pH from acidic to basic were tested. The results showed that varying ranges of pH did not affect the test results.

### E. Interference












The compounds listed below showed no interference.

Acetaminophen	Doxepin	Nortriptyline
Acetone (1000mg/dL)	Ecgonine methyl ester	Noscapine

Acetophenetidin	Ephedrine	O-Hydroxyhippuric acid
Acetylsalicylic acid	Erythromycin	Octopamine
Albumin (100mg/dL)	Ethanol (1%)	Oxalic acid (100 mg/dL)
Albuterol	Fenoprofen	Oxazepam
Aminopyrine	Fluphenazine	Oxolinic acid
Amitriptyline	Furosemide	Oxymetazoline
Amobarbital	Galactose (10mg/dL)	Papaverine
Amoxicillin	Gamma Globulin (500mg/dL)	Penicillin G
Ampicillin	Gentisic acid	Perphenazine
Apomorphine	Glucose (3000mg/dL)	Phencyclidine
Ascorbic acid	Hemoglobin	Phenelzine
Aspartame	Hydralazine	Phenobarbital
Atropine	Hydrochlorothiazide	Prednisone
Benzilic acid	Hydrocortisone	Propoxyphene (50ug/ml)
Benzoic acid	Hydroxytyramine	Propranolol
Benzoyllecgonine	Ibuprofen	Pseudoephedrine
Bilirubin	Imipramine	Quinine
Boric Acid (1%)	Isoproterenol	Ranitidine
Bupropion	Isoxsuprine	Riboflavin (10mg/dL)
Caffeine	Ketamine	Salicylic acid
Carbamazepine	Ketoprofen	Secobarbital
Chloral hydrate	Labetalol	Serotonin (5-Hydroxytyramine)
Chloramphenicol	Lidocaine	Sulfamethazine
Chlorothiazide	Loperamide	Sulindac
Chlorpromazine	Maprotiline	Tetrahydrocortisone3-(β-Dglucuronide)
Cholesterol	Meperidine	Tetrahydrocortisone 3-acetate
Clomipramine	Meprobamate	Tetrahydrozoline
Clonidine	Methapyrilene	Thiamine
Cortisone	Methaqualone	Thioridazine
Cotinine	Methoxyphenamine	Triamterene
Creatinine	Metronidazole (300ug/ml)	Trifluoperazine
Cyclobenzaprine	N-Acetylprocainamide	Trimethoprim
Deoxycorticosterone	NaCl (4000mg/dL)	Tyramine
Desipramine	Nalidixic acid	Urea (2000mg/dL)

Dextromethorphan	Naloxone	Uric acid
Diclofenac	Naltrexone	Valproic acid (250ug/ml)
Diflunisal	Naproxen	Venlafaxine
Digoxin	Niacinamide	Verapamil
Diphenhydramine	Nicotine	Zomepirac
DL-Tryptophan	Nifedipine	β-Estradiol
DL-Tyrosine	Norethindrone	

**Index of Symbols**

	Consult instructions for use or consult electronic instructions for use		Contains sufficient for <n> tests		Store between 35.6-86°F (2-30°C)
	<i>In vitro</i> diagnostic medical device		Batch code		Catalogue number
	Manufacturer		Use-by date		Do not re-use
	Do not use if package is damaged and consult instructions for use		Unique device identifier		

**CLIA WAIVED**

Manufactured for:  
**EKLA Corporation**  
1707 Quincy Ave. Suite 127, Naperville, IL 60540 USA  
Customer Service Phone: 1-800-320-4215/1-630-258-6242