

novaplus[™] hCG Rapid Pregnancy Test



- Easy to use, minimal training required.
- One-step, no additional reagents required.
- Proven accuracy of >99%.
- Compact design maximizes storage space in laboratory and Point Of Care locations.
- Includes hCG variance H-hCG and hCGBcf, for the detection of early and mid-stage pregnancy.
- Combo Test, cut off 10 mIU/hCG serum, 20 mIU/hCG urine, can detect 1×10^6 6 pmol/ml level of hCGBcf and/or H-hCG.



Item No.	Description	Analytical Sensitivity	Eaches/case	CPT Code
VPT1-CASS	Urine hCG Cassette with Pipette	20 mIU/hCG Urine	40	81025
VPT2-CASS	Combo Urine/Serum hCG Cassette with Pipette	10 mIU/hCG Serum 20 mIU/hCG Urine	40	84703

hCG Rapid Pregnancy Test for the determination of hCG in urine and serum specimens. This test kit is used to obtain a visual qualitative result for early detection of pregnancy.

